



Internal Training Only



This device is breast massager

BELUNA will help you to maintain the healthy breast

 <p>Ergonomic 3D Design</p>	 <p>Treatment device with powerful motor</p>	 <p>1:1 Treatment for the desired effect</p>	 <p>Solely product in South Korea with three circulation modes</p>
--	---	---	---



TOTAL SOLUTION FOR HEALTHY BREAST

BELUNA

Hold ON!

Are you confused if this product is for you or not?
Check this out first!

1. AFTER/BEFORE SURGERY

- For those who need to ease the breast skin tension before surgery
- For those who want to stimulate the image after surgery
- For those who want a quick recovery and natural breast line after surgery

2. BREAST VOLUME UP

- For those who need volume up breast on special occasion
- For those who want a perfect smooth fit and S-Line
- For those who want a natural and comfortable volume up by increasing the attach feeling of breast on the skin

3. ELASTICITY IMPROVEMENT/LINE MAINTENANCE

- For those who concern about saggy breast after dieting
- For those who need confidence boost after breast feeding
- For those who concern about breast fat, accessory breast

Does one of them apply to you? Schedule your consultation NOW!

CG Bio Co., Ltd.
3rd FL, 211, Itaewon-ro, Yongsan-gu, Seoul, Korea
T. 02-550-8300 | www.cgbio.co.kr

CGBIO
CELL & GROWTH FACTOR
BIOTECHNOLOGY



1. How BELUNA works?

BELUNA is using a negative pressure to stimulate breast smoothly



- ① **Manage the flab around breast**
Gather the breast fat, accessory breast around breast together
- ② **Put the strained breast at ease**
Put the strained and suffocated breast that caused by bra at ease
- ③ **Fill it UP, volume it UP**
Possibility to use the device according to your needs by using negative pressure

2. How to use BELUNA?

You can choose any BELUNA mode that applied to you when using it.
You can adjust the device pressure based on your condition

Continuous Mode

This mode constantly provides negative pressure while using the device.



PURPOSE OF USAGE
after breast augmentation

USAGE TIME
min. 1 hour of use each day

USAGE PERIOD
min. 2 weeks of use

Cyclic Mode

This mode sets the lowest and highest pressures to maintain the two pressures alternately for a certain period of time and operate in circulation.



PURPOSE OF USAGE
Temporary breast enlargement before breast augmentation (fat transplantation/implant transplantation)

USAGE TIME
min. 1 hour of use each day

USAGE PERIOD
min. 2 weeks of use

Massage Mode

This mode operates by setting two pressures, either the lowest or highest, without restricted time.



PURPOSE OF USAGE
Improving the touch feeling, preventing pain, and tightening the breast.

USAGE TIME
min. 1 hour of use each day

Usage period
min. 2 weeks of use

3. Why BELUNA? (Q&A)

Q. Why do we need BELUNA before implant transplantation?

A. When implants are inserted inside the breast without sufficient space, it can result in high risk of side effects including the visibility of implant edge (rippling), capsular contracture, etc. Using BELUNA can help to reduce the likelihood these kinds of side effects from surgery by creating enough space for implant in the breast tissue.

Q. Why do we need BELUNA after implant transplantation?

A. Using BELUNA after surgery can help shaping the implant transplantation result into more natural look and improvement of touch-feeling through stimulation of massage.

Q. Why do we need BELUNA before fat grafting?

A. Growth rate of the fat grafting is the most important! Enough space must be secured to place the fat so the attachment rate also increases. Using BELUNA before fat grafting can increase the tissue to secure enough space to place the fat with stimulating vascular production and blood flow to provide enough nutrients for grafted fat cells.

Q. Why do we need BELUNA after fat grafting?

A. Using BELUNA after fat grafting can accelerate waste discharge through lymph nodes. Also, accelerating vascular production and activating blood circulation might help to increase the growth rate of grafted fat by facilitating the supply of nutrients and hormones.

FACT 1

"Have you ever seen a picture of Native Africans who wear many neck rings to fill their necks? They add necklace one by one every certain time to stretch their neck. This phenomenon shows that our body has properties that can gradually transform according to continuous stimulation that has applied to it. Same thing also applied to person with a short stature or people who need a correction of different legs length due to accidents. By applying the Illizarov technique, their leg bones are intentionally cut and widened 1mm every day so the leg tissues (e.g., bones and blood vessels) will grow up to 10 after a year. This technique can be seen in the same way as BELUNA works."



FACT 2

BELUNA uses the negative pressure by applying NPWT system.

What is NPWT (Negative Pressure Wound Therapy)? Negative pressure collects the blood flow needed for cell proliferation that will stimulates cells around the wound to induce the production of collagen which has been lost as the skin ages. BELUNA is a product developed and designed based on CGBIO's Curasys technology. NPWT (Negative Pressure Wound Therapy) has already proven in professional literature for their effectiveness.

